

## Cement Removal

QwikStrip™ Serrated Strips are designed to make crown and bridge clean-up quick, safe, and simple after cementation of crowns, inlays, onlays, bridges or veneers. Cleaning up is never this easy!

\* Hold the Serrated Strip with your thumb and index finger and insert the serrated blade between the teeth with a gentle "sawing motion" to easily remove the overflow of cement.

## Finishing, Polishing, Stain Removal or Overhang Correction

QwikStrip™ single-sided (straight or curved) version can be used in the anterior or the posterior of the mouth. Finish or polish any restoration to perfection with the QS-SF (superfine) strips. Using the curved version (QSC-SF) will enable you to remove stains, preserve the anatomical contour where needed, or to correct overhangs.

## Interproximal Reduction (IPR)

When used for any IPR procedure, simply follow the QwikStrip™ color coded sequence guide provided for achieving precise Interproximal reduction (IPR). The guide will enable you to achieve from .07 mm up to .5mm of space. Use all versions of QwikStrip™ (single-sided, double-sided, and curved) in any area of the mouth.

\* Following the sequence guide will ensure that you accomplish the desired result. For safety, use floss through the tether holes on the QwikStrips.

## IPR Sequence Guide

In order to create precise interproximal spaces for any Orthodontic procedure (including Invisalign™, Clear Correct™, Six Month Smiles™, or MTM™) use **QwikStrip™** color coded strips according to the sequence guide outlined below. The reduction process (IPR) should begin with the yellow single-sided grit (.07 mm), and continue in sequence until the desired space is achieved up to .5mm.

**Check the reduction with an approved measurement guide before moving to the next color (grit) in the sequence.**

*Note: Use all versions of QwikStrip™ (single-sided, double-sided, or curved) in any area of the mouth. However, for maintaining the anatomical contour on posterior teeth or when refining the natural anatomy of the anterior teeth, finish by using the curved version.*

All widths are ± .05 mm.

