

Oblouky pro Alexander LTS

- Vhodné při použití zámků v preskripci Alexander LTS
- Ze všech nabízených tvarů LTS oblouky mají nejužší vzdálenost mezi špičkami

Alexander LTS Wire Recommendations

		Maxillary		Mandibular		Treatment Phase
.017 x .025	Therma-Ti Lite	857-943A	W	857-953A	W	Initial treatment wire
.016	Force I NiTi	857-506A	W	857-516A	W	Initial treatment wire
.016 x .022	8 Stranded	857-902A		857-912A		OPTIONAL Initial treatment wire
.017 x .025	Twist Wire	857-236A		857-256A		OPTIONAL Initial treatment wire
.016	Stainless Steel	853-236A	W	853-246A	W	Transitional non-extract space closure
.016 x .022	Stainless Steel	853-201A	W	853-221A	W	Transitional arch leveling
.017 x .025	Force I NiTi	857-537A	W	857-547A	W	Transitional arch leveling
.017 x .025	Beta Titanium	857-672A	W	857-682A	W	Transitional arch leveling
.019 x .025	Beta Titanium	857-674A	W	857-684A	W	Transitional arch leveling - .022 user
.017 x .025	Stainless Steel	853-204A	W	853-224A	W	Transitional Extraction Space Closure/Finishing
.019 x .025	Stainless Steel	853-207A	W	853-227A	W	Transitional Extraction Space Closure/Finishing - .022 User

