



Innovation. Quality. Efficiency.

PRODUCTS CONTAINING ANTIMICROBIAL ADDITIVE

03/27/2014

Zirc incorporates antimicrobial compounds into the following products: Tubs & accessories, Trays, Instrument Holders, Bur/Endo Holders, Mouth Mirrors, Implant and Restorative Organizers, E-Z Access Shelf and Accessories, Lab Pans, Evacuation Screens, Steri Soaker, Water Wise, Insti Dam Dispensor, Waste Drop Rings, and Handi Hopper.

Antimicrobial compounds are intended to protect parts against staining, odors, and degradation due to fungi, mold, yeasts, mildew, algae and other microbes. The incorporation of antimicrobial compounds in our plastic products is not intended to protect users or others against food-borne or disease-causing bacteria. For more information please reference EPA Registration Notice 2000-1.

Zirc Dental Products, Inc.

3918 Highway 55 S.E. • Buffalo, MN 55313
Tel. 763-682-6636 ~ 800-328-3899 ~ Fax: 763-682-6604
www.zirc.com